



is now part of **ZIMMER BIOMET**
 Your progress. Our promise.®

EUROPE

Parc d'entreprises du Grand Troyes

5 rue de Berlin

10300 Sainte-Savine, France

Adresse postale : CS 80002

10302 Sainte-Savine CEDEX

+33 (0)3 25 82 32 63

ASIA-PACIFIC

Suite 1001, Great Eagle Centre,

23 Harbour Road,

Wanchai,

Hong Kong

+ (852) 3955 1838

LATIN AMERICA

Av. Pereira Barreto, 1395

19º andar

Torre Sul - Bairro Paraiso

Santo André - São Paulo

CEP: 09190-610

Brazil

+ 55 11 43327755

UNITED STATES

13785 Research Boulevard

Suite 200

Austin, TX 78750

USA

+1 (512) 344 3333

www.ldr.com



MC+®
 CERVICAL MODULAR CAGE



Document intended for the exclusive use of healthcare professionals. MC+® is a range of CE marked medical devices of classes IIb, IIa and I manufactured by the LDR Medical company. The conformity assessment of classes IIb and IIa devices was carried out by the Notified Body LNE/G-MED N°0459. The MC+® implant is a sterile cage intended for an arthrodesis of the cervical vertebrae through anterior approach. The MC+® instruments constitute an instrument set intended to allow the MC+® implant implantation. Before any surgical procedure, read carefully the instructions and the surgical technique.

©2016 LDR Holding Corporation and its affiliates. All Rights Reserved. Corporation or its affiliates in France, the United States or other countries.

LDR, LDR Spine, LDR Medical, Avenue, Bt-Pack, Bf+, Bf+(ph), C-Plate, EasySpine, InterBRIDGE, MC+, MIVO, Mobl, Mobl-C, Moblisc, ROI, ROI-A, ROI-C, ROI-T, SpineTune, and VerteBRIDGE are trademarks or registered trademarks of LDR Holding Corporation.



REF : MC B 1 EN 05.2017 D

Non contractual photos. LDR reserves the right to modify its products without notice.

MC+[®]

CERVICAL MODULAR CAGE

A design supporting stability

- The anatomical shape of the MC+[®] cage (upper convex endplate in sagittal and frontal planes) allows restoration of disc height as well as intersomatic space lordosis while offering excellent primary stability to the implant.
- Both the superior and inferior surfaces of the cage feature retaining ridges limiting the risk of implant migration.
- The large weight bearing surface peripheral to the bone graft area increases primary stability and allows for immediate loading.
- The self-retaining anchoring clip supplements the primary stability of the cage by inserting through the superior endplate of the lower vertebrae.
- The exclusive design of the anchoring clip provides its self-locking within the cage, easy impaction, and reliable supplemental fixation without weakening the vertebral endplate.



Patient safety optimized

- The MC+[®] cage is in PEEK-Optima[®] with a Young's modulus close to that of bone. Its radiolucency and the use of tantalum markers allow post-operative visualization and verification of the cage positioning and fusion.
- The intersomatic cages MC+[®] are implanted using a reproducible surgical technique.
- The MC+[®] cage offers a wide unique graft surface with the possibility to use either an anatomical bone substitute or autologous bone graft to perform the fusion.



A wide range that fits to patients anatomy

The large range of cages (footprint, height) allows adaptation of the implant to varying patient anatomy, and makes multi-level surgeries supportable.

Product range and sizing *

Footprints mm (depth x width)	12x14
	12x15,5
	14x14
	14x15,5
	14x17
Heights	4,5mm ^{**} ; 5mm; 6mm; 7mm; 8mm
Anchoring clips	One-size-fits-all

* Product availability may vary by country and market.
** Height not available for Footprint 14x17.



- The insertion of the anchoring clip is self-guided through the cage holder in the angular direction (50° angle in regards to the axis of the disc) as well as in depth allowing for perfect positioning (impactor with stop).

- All MC+[®] implants are delivered sterile with sterility indicator.

Compact instrumentation

The MC+[®] instruments are simple and functional and have been designed by and for spinal surgeons.

It allows for a simple, optimized and reproducible surgical technique.

